



The Challenge

With packaging now a core element in the marketing mix, pharmaceutical companies seek packaging solutions that provide a positive patient experience, are child resistant yet user friendly, and support medication adherence and brand loyalty.

In addition, with the vast majority of drug formulations requiring protection from moisture and/or oxygen to meet stability requirements, the packaging might require passive as well as active barrier solutions which help to extend shelf life and ensure regulatory compliance.



The Solution ActiveGuard™

LOG Pharma Packaging's ActiveGuard™ meets all of these challenges. The patient-centric solution is prefilled with moisture desiccants and/or oxygen absorbers according to customer requirements, and placed firmly in the neck of LOG's specially designed barrier bottles.

Features & Benefits:

Patient Centricity	Health Improvement	Packaging Optimization
<ul style="list-style-type: none"> ✔ Eliminates the need for desiccant sachets or canisters ✔ Cannot fall in or out of the packaging ✔ Removes the risk of swallowing the sachet/canister accidentally ✔ Allows controlled dispensing of tablets 	<ul style="list-style-type: none"> ✔ Allows the use of heat seal liners for hermetic sealing and tamper evidence ✔ Amount of desiccants and/or oxygen absorbers according to customer requirements ✔ No direct contact with tablets, avoiding possible damage 	<ul style="list-style-type: none"> ✔ Positioned in previously unused space ✔ Does not reduce valuable tablet-space ✔ Reduces head space ✔ Easy positioning validation

Neck finishes: Available for **33mm** (fillable with up to 2gr desiccant) and **38mm** (fillable with up to 4gr desiccant)

ActiveGuard™

IT'S ALL ABOUT THE

Click

The most user friendly packaging solution for moisture sensitive formulations



For more information & samples, please contact:

LOG Pharma Packaging
T. +972 4 6756338 | sales@logpac.com | www.logpac.com